


**Product Name:** Neopentylboronic acid

**Catalog Number:** N1152



**Sizes Available:** 1 g, 5 g, and larger sizes available

**Molecular weight:** 115.97 g/mol

**Molecular Formula:** C<sub>5</sub>H<sub>13</sub>BO<sub>2</sub>

**CAS Number:** 701261-35-0

**Storage:** Store at -20 °C, under dry conditions.

**Synonyms:** NA

**Uses:** Transition Metal Coupling Reactions

## Selected References

### Boronic Acid Derivative Coupling Reactions

- Altenhoff, Gereon, Richard Goddard, Christian W. Lehmann, and Frank. Glorius. Sterically Demanding, Bioxazoline-Derived N-Heterocyclic Carbene Ligands with Restricted Flexibility for Catalysis. *Journal of the American Chemical Society* 126, no. (2004): 15195–201. <https://doi.org/10.1021/ja045349r>.
- Arvela, Riina K., Nicholas E. Leadbeater, Michael S. Sangi, Victoria A. Williams, Patricia Granados, and Robert D. Singer. A Reassessment of the Transition-Metal Free Suzuki-Type Coupling Methodology. *Journal of Organic Chemistry* 70, no. (2005): 161–68. <https://doi.org/10.1021/jo048531j>.
- Billingsley, Kelvin, and Stephen L. Buchwald. Highly Efficient Monophosphine-Based Catalyst for the Palladium-Catalyzed Suzuki-Miyaura Reaction of Heteroaryl Halides and Heteroaryl Boronic Acids and Esters. *Journal of the American Chemical Society* 129, no. (2007): 3358–66. <https://doi.org/10.1021/ja068577p>.
- Billingsley, Kelvin L., Kevin W. Anderson, and Stephen L. Buchwald. A Highly Active Catalyst for Suzuki-Miyaura Cross-Coupling Reactions of Heteroaryl Compounds. *Angewandte Chemie, International Edition* 45, no. (2006): 3484–88. <https://doi.org/10.1002/anie.200600493>.
- Gillis, Eric P., and Martin D. Burke. A Simple and Modular Strategy for Small Molecule Synthesis: Iterative Suzuki-Miyaura Coupling of B-Protected Haloboronic Acid Building Blocks. *Journal of the American Chemical Society* 129, no. (2007): 6716–17. <https://doi.org/10.1021/ja0716204>.
- Herradura, Prudencio S., Kathleen A. Pendola, and R. Kiplin. Guy. Copper-Mediated Cross-Coupling of Aryl Boronic Acids and Alkyl Thiols. *Organic Letters* 2, no. (2000): 2019–22. <https://doi.org/10.1021/ol005832g>.
- Kataoka, Noriyasu, Quinetta Shelby, James P. Stambuli, and John F. Hartwig. Air Stable, Sterically Hindered Ferrocenyl Dialkylphosphines for Palladium-Catalyzed C-C, C-N, and C-O Bond-Forming Cross-Couplings. *Journal of Organic Chemistry* 67, no. (2002): 5553–66. <https://doi.org/10.1021/jo025732j>.
- Kinzel, Tom, Yong Zhang, and Stephen L. Buchwald. A New Palladium Precatalyst Allows for the Fast Suzuki-Miyaura Coupling Reactions of Unstable Polyfluorophenyl and 2-Heteroaryl Boronic Acids. *Journal of the American Chemical Society* 132, no. (2010): 14073–75. <https://doi.org/10.1021/ja1073799>.
- Knapp, David M., Eric P. Gillis, and Martin D. Burke. A General Solution for Unstable Boronic Acids: Slow-Release Cross-Coupling from Air-Stable MIDA Boronates. *Journal of the American Chemical Society* 131, no. (2009): 6961–63. <https://doi.org/10.1021/ja901416p>.
- Liebeskind, Lanny S., and Jiri. Srogl. Thiol Ester-Boronic Acid Coupling. A Mechanistically Unprecedented and General Ketone Synthesis. *Journal of the American Chemical Society* 122, no. (2000): 11260–61. <https://doi.org/10.1021/ja005613q>.
- McGuinness, David S., and Kingsley J. Cavell. Donor-Functionalized Heterocyclic Carbene Complexes of Palladium(II): Efficient Catalysts for C-C Coupling Reactions. *Organometallics* 19, no. (2000): 741–48. <https://doi.org/10.1021/om990776c>.
- Miyaura, N., T. Yanagi, and A. Suzuki. The Palladium-Catalyzed Cross-Coupling Reaction of Phenylboronic Acid with Haloarenes in the Presence of Bases. *Synthetic Communications* 11, no. (1981): 513–19. <https://doi.org/10.1080/00397918108063618>.

Miyaura, Norio. Organoboron Compounds. *Topics in Current Chemistry* 219, no. (2002): 11–59.

Nguyen, Hanh Nho, Xiaohua Huang, and Stephen L. Buchwald. The First General Palladium Catalyst for the Suzuki-Miyaura and Carbonyl Enolate Coupling of Aryl Arenesulfonates. *Journal of the American Chemical Society* 125, no. (2003): 11818–19. <https://doi.org/10.1021/ja036947t>.

Qiao, Jennifer X., and Patrick Y. S. Lam. Copper-Promoted Carbon-Heteroatom Bond Cross-Coupling with Boronic Acids and Derivatives. *Synthesis*, no. (2011): 829–56. <https://doi.org/10.1055/s-0030-1258379>.

Seiple, Ian B., Shun Su, Rodrigo A. Rodriguez, Ryan Gianatassio, Yuta Fujiwara, Adam L. Sobel, and Phil S. Baran. Direct C-H Arylation of Electron-Deficient Heterocycles with Arylboronic Acids. *Journal of the American Chemical Society* 132, no. (2010): 13194–96. <https://doi.org/10.1021/ja1066459>.

Shirota, Yasuhiko, Motoi Kinoshita, Tetsuya Noda, Kenji Okumoto, and Takahiro. Ohara. A Novel Class of Emitting Amorphous Molecular Materials as Bipolar Radical Formants: 2-{4-[Bis(4-Methylphenyl)Amino]Phenyl}-5-(Dimesitylboryl)Thiophene and 2-{4-[Bis(9,9-Dimethylfluorenyl)Amino]Phenyl}-5-(Dimesitylboryl)Thiophene. *Journal of the American Chemical Society* 122, no. (2000): 11021–22. <https://doi.org/10.1021/ja0023332>.

Stambuli, James P., Ryoichi Kuwano, and John F. Hartwig. Unparalleled Rates for the Activation of Aryl Chlorides and Bromides: Coupling with Amines and Boronic Acids in Minutes at Room Temperature. *Angewandte Chemie, International Edition* 41, no. (2002): 4746–48. <https://doi.org/10.1002/anie.200290036>.

Tsuboyama, Akira, Hironobu Iwawaki, Manabu Furugori, Taihei Mukaide, Jun Kamatani, Satoshi Igawa, Takashi Moriyama, et al. Homoleptic Cyclometalated Iridium Complexes with Highly Efficient Red Phosphorescence and Application to Organic Light-Emitting Diode. *Journal of the American Chemical Society* 125, no. (2003): 12971–79. <https://doi.org/10.1021/ja034732d>.

Wang, Dong-Hui, Masayuki Wasa, Ramesh Giri, and Jin-Quan. Yu. Pd(II)-Catalyzed Cross-Coupling of Sp<sup>3</sup> C-H Bonds with Sp<sup>2</sup> and Sp<sup>3</sup> Boronic Acids Using Air as the Oxidant. *Journal of the American Chemical Society* 130, no. (2008): 7190–91. <https://doi.org/10.1021/ja801355s>.

Wright, Stephen W., David L. Hageman, and Lester D. McClure. Fluoride-Mediated Boronic Acid Coupling Reactions. *Journal of Organic Chemistry* 59, no. (1994): 6095–97. <https://doi.org/10.1021/jo00099a049>.

Hazardous Properties and Cautions: The toxicological and pharmacological properties of this compound are not fully known. For further information see the MSDS on request. **Neopentylboronic acid** is manufactured, shipped according to standard practices, and intended for research and development in a laboratory utilizing prudent procedures for handling chemicals of unknown toxicity, under the supervision of persons technically qualified to evaluate potential risks and authorized to enforce appropriate health and safety measures. As with all research chemicals, precautions should be taken to avoid unnecessary exposures or risks.

Warranty and Disclaimer: Frontier Scientific, Inc. warrants the product conforms to the specifications stated herein. In the event of nonconformity, Frontier will replace products or refund purchase price, at its sole option, and Frontier shall not be responsible for any other loss or damage, whether known or foreseeable to Frontier. No other warranties apply, express or implied, including but not limited to warranty of fitness for any purpose or implied warranty of merchantability. Purchaser is solely responsible for all consequences of its use of the product and Frontier assumes no responsibility therefore, ) including success of purchaser's research and development, or health or safety of any uses of the product.

