


Product Name: 6-Bromoimidazo[1,5-a]pyridine

Catalog Number: Y15132



Sizes Available: 1 g, 5 g, 25 g, and larger sizes available

Molecular weight: 197.03 g/mol

Molecular Formula: C₇H₅BrN₂

CAS Number: 1239880-00-2

Storage: Store at 2-8 C°, under dry conditions.

Synonyms: 6-bromoimidazo[1,5-a]pyridine, imidazo[1,5-a]pyridine, 6-bromo-

Uses: Synthesis building block, Organic Synthesis, imidazole-pyridine nitrogen heterocycle, synthesis, bromine reactions

6-Bromoimidazo[1,5-a]pyridine, is a synthetic fine chemical useful in the synthesis of pharmaceuticals and fine organic chemicals.

Selected References:

1. Albaugh, P. A.; Choi, H.-S.; Chopiuk, G.; Fan, Y.; Rucker, P. V.; Wang, Z. Preparation of imidazopyridazines and similar compounds as therapeutic kinase inhibitors. [WO2009140128A2, 2009.](#)
2. Almirante, L.; Mugnaini, A.; De Toma, N.; Gamba, A.; Murmann, W.; Hidalgo, J., [Imidazole derivatives. IV. Synthesis and pharmacologic activity of oxygenated derivatives of imidazo\[1,2-a\]pyridine.](#) *J. Med. Chem.* **1970**, 13 (6), 1048-51.
3. Bruice, T. C.; Schmir, G. L., [Imidazole catalysis. II. The reaction of substituted imidazoles with phenyl acetates in aqueous solution.](#) *J. Am. Chem. Soc.* **1958**, 80, 148-56.
4. Gudmundsson, K.; Boggs, S. D. Preparation of imidazo[1,2-a]pyridine derivatives as chemokine receptor ligands with anti-HIV activity. [WO2006026703A2, 2006.](#)
5. Ismail, M. A.; Arafa, R. K.; Wenzler, T.; Brun, R.; Tanious, F. A.; Wilson, W. D.; Boykin, D. W., [Synthesis and antiprotozoal activity of novel bis-benzamidino imidazo\[1,2-a\]pyridines and 5,6,7,8-tetrahydro-imidazo\[1,2-a\]pyridines.](#) *Bioorg. Med. Chem.* **2008**, 16 (2), 683-691.
6. Ni, Z.-J.; Pecchi, S.; Burger, M.; Han, W.; Smith, A.; Atallah, G.; Bartulis, S.; Frazier, K.; Verhagen, J.; Zhang, Y.; Iwanowicz, E.; Hendrickson, T.; Knapp, M.; Merritt, H.; Voliva, C.; Wiesmann, M.; Legrand, D. M.; Bruce, I.; Dale, J.; Lan, J.; Levine, B.; Costales, A.; Liu, J.; Pick, T.; Menezes, D. Preparation of imidazo[1,2-a]pyridines and imidazo[1,2-b]pyridazines as PI-3 kinase inhibitors. [WO2007095588A1, 2007.](#)
7. Ratcliffe, A. J.; Alam, M.; Beever, R. E.; Davenport, R. J.; Davies, N.; Haughan, A. F.; Jones, M. W.; Lowe, C.; Perry, B. G.; Phillips, D. J.; Pitt, W. R.; Sharpe, A. Preparation of aminopyrimidines as JNK inhibitors. [WO2006038001A1, 2006.](#)
8. Ratni, H.; Karp, G. M.; Weetall, M.; Naryshkin, N. A.; Paushkin, S. V.; Chen, K. S.; McCarthy, K. D.; Qi, H.; Turpoff, A.; Woll, M. G.; Zhang, X.; Zhang, N.; Yang, T.; Dakka, A.; Vazirani, P.; Zhao, X.; Pinard, E.; Green, L.; David-Pierson, P.; Tuerck, D.; Poirier, A.; Muster, W.; Kirchner, S.; Mueller, L.; Gerlach, I.; Metzger, F., [Specific Correction of Alternative Survival Motor Neuron 2 Splicing by Small Molecules: Discovery of a Potential Novel Medicine To Treat Spinal Muscular Atrophy.](#) *J. Med. Chem.* **2016**, 59 (13), 6086-6100.
9. Uchikawa, O.; Sakai, N.; Terao, Y.; Suzuki, H. Preparation of fused heterocyclic compounds as apoptosis signal regulating kinase 1 (ASK1) inhibitors. [WO2008016131A1, 2008.](#)

10. Zhou, X.; Yan, H.; Ma, C.; He, Y.; Li, Y.; Cao, J.; Yan, R.; Huang, G., [Copper-Mediated Aerobic Oxidative Synthesis of 3-Bromo-imidazo\[1,2-a\]pyridines with Pyridines and Enamides](#). *J. Org. Chem.* **2016**, *81* (1), 25-31.

Hazardous Properties and Cautions: The toxicological and pharmacological properties of this compound are not fully known. For further information see the MSDS on request. **6-Bromoimidazo[1,5-a]pyridine** is manufactured, shipped according to standard practices, and intended for research and development in a laboratory utilizing prudent procedures for handling chemicals of unknown toxicity, under the supervision of persons technically qualified to evaluate potential risks and authorized to enforce appropriate health and safety measures. As with all research chemicals, precautions should be taken to avoid unnecessary exposures or risks.

Warranty and Disclaimer: Frontier Scientific, Inc. warrants the product conforms to the specifications stated herein. In the event of nonconformity, Frontier will replace products or refund purchase price, at its sole option, and Frontier shall not be responsible for any other loss or damage, whether known or foreseeable to Frontier. No other warranties apply, express or implied, including but not limited to warranty of fitness for any purpose or implied warranty of merchantability. Purchaser is solely responsible for all consequences of its use of the product and Frontier assumes no responsibility therefore,) including success of purchaser's research and development, or health or safety of any uses of the product.